

Now Get More With Less...

Seven-fold Reduction in Sample Cell Volume (200 μl)

Introducing the new iTC₂₀₀ from MicroCal

An evolution in ITC instrumentation, the iTC₂₀₀ is more sensitive, faster, and easier to use than the VP-ITC, which has been the industry standard.



Advantages of the iTC₂₀₀ vs. the VP-ITC:

- Sample Usage: The iTC₂₀₀ requires up to 7x less sample material than the VP-ITC.
- Volume Usage: The iTC₂₀₀ requires 7x less sample volume than the VP-ITC while providing the same great results. The sample cell volume of the iTC₂₀₀ is only 200 μl and the syringe volume is 40 μl .
- Throughput: With faster equilibration times, the iTC₂₀₀ is two to four times faster than the VP-ITC. The iTC₂₀₀ can also be upgraded to provide sample throughput of up to 50 samples per day with a capacity to process as many as 384 samples unattended.
- Improved sample loading and cell cleaning.
- The iTC₂₀₀ is controlled by a new intelligent user interface that facilitates experimental design and automatically processes data at the end of sample runs. Data is presented in Excel format for further analysis or data transfer.

While sample and volume usage are reduced, the iTC₂₀₀ still provides the same great data you have come to expect with the VP-ITC.

As can be seen in Table 1, the results for K and ΔH are in good agreement between the VP-ITC and iTC₂₀₀. However, the iTC₂₀₀ requires up to 7x less protein.

Table 1. RNase Quantity Requirements Based on c value Requirements

Values	VP-ITC (1,400 μ L cell volume)				iTC ₂₀₀ (200 μ L cell volume)			
c value	50	25	10	5	50	25	10	5
Protein Concentration (mM)	0.035	0.018	0.0070	0.0035	0.035	0.018	0.0070	0.0035
μ g of Protein	685	343	137	69	98	49	20	10
Experimental Results								
$K_A \times 10^6$ (M ⁻¹)	1.5	1.5	1.6	1.6	1.6	1.6	1.6	1.9
K_D (μ M)	0.66	0.66	0.62	0.60	0.62	0.59	0.61	0.52
ΔH (kcal/mol)	-18	-17	-16	-16	-17	-17	-15	-15

Quantities calculated for the ITC titration of RNase-2'CMP assuming $K_A = 1.4 \times 10^6$, $\Delta H = -15$ kcal/mol, MW = 13,700. Twelve injections were done per experiment.

The iTC₂₀₀ is much faster than the VP-ITC with faster equilibration times.

As can be seen in Table 2, two runs per hour can be accomplished with the iTC₂₀₀ compared to the VP-ITC which takes over an hour to complete a single run. Sample rates are doubled with the iTC₂₀₀ and can be doubled again by reducing the number of injections from the traditional 24 to 12. By implementing a reduced injection strategy, up to four runs per hour can be accomplished with the iTC₂₀₀.

Table 2. ITC Sample Processing Time

Function	VP-ITC		iTC ₂₀₀	
# of injections	24	12	24	12
Cleaning (minutes)	8	8	8	8
Loading (minutes)	2	2	2	2
Equilibration (minutes)	20	20	4	4
Titration (minutes)	96	48	36	18
TOTAL (minutes)	126	78	50	32

Ultrasensitive Calorimetry for the Life Sciences™

MicroCal, LLC
 22 Industrial Drive East, Northampton, MA 01060
 United States of America
 T +1 413-586-7720 / 1-800-633-3115
 F +1 413-586-0149
 E info@microcal.com
 www.microcal.com

MicroCal, LLC Europe
 2 Warren Yard, Warren Farm Office Village
 Wolverton Mill, Milton Keynes MK12 5NW, United Kingdom
 T +44 1908 576330
 F +44 1908 576339
 E info@microcal.eu.com
 www.microcal.com

